SMARTRING

QUARK FitRing

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[clic aquí](#) | [click here](#)

ES
ESPAÑOL

EN
ENGLISH

www.foxbox.com.ar

MANUAL DE USUARIO FITRING

Uso apropiado

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USER'S MANUAL FITRING

Proper use

Please wear the ring on the appropriate finger, choosing a comfortable size that is neither loose nor tight. It is recommended to wear it on the index, middle and/or ring finger of non-dominant hands to better control the data. When taking any measurements, the sensors inside the ring should be firmly attached to the finger.

Download the App

Scan the QR code to download and install the app. The smart ring must be linked to the app before use.

Click here (or Scan the QR):

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

DISPONIBLE EN
Google Play

Product description


Charging case


Magnetic Charger

Scheme of use


Device connection

1. Charge the ring before the first use.
2. Download the Da Rings application on your phone from the official store (App Store or Google Play).
3. Activate Bluetooth on your phone.
4. Open the Da Rings app and fill in the requested information to set up your profile.
5. Add the "FOXBOX | FRG" ring and pair the device with the application. Make sure that Bluetooth is turned on before this step.
6. Put on the ring to start measuring and recording data automatically.

Charging and ignition

Connect the ring with the magnetic charging cable, the red light turns on and the green light flashes to turn on.

Note: The red light flashes at intervals during charging. The green light flashes when fully charged.

Ring functions

Health records

The ring monitors in real time various health indicators, including:

Sleep: analysis of sleep phases and quality of rest.

Heart rate: Continuous heart rate measurement.

Blood oxygen: Oxygen saturation levels.

Stress level: Assessment based on heart rate variability.

Menstrual cycle tracking: Cycle phase recording and predictions.

Activity logs (18 modes)

The ring allows you to record various sports and physical activities:

Aerobic and endurance: walking, running, cycling, stationary cycling.

High intensity: Skipping, fitness, free training.

Sports: Badminton, basketball, soccer, climbing, tennis, rugby, golf, yoga.

Other: Dancing, elliptical, rowing.

GPS Activity Logging (5 modes)

The ring uses the phone's GPS to improve accuracy in certain activities:

GPS walking

GPS running

GPS cycling

GPS mountain running

GPS steps

Basic Specifications

Memory Storage: **1Mb ROM + 256Kb RAM**

Bluetooth version: **5.1**

Case Battery: **220 mAh**

Battery Life: **4 Days**

Charging Time: **2.5 Hs.**

Water resistance: **5ATM**

Precautions

Before using and operating the equipment, read and follow the following precautions to ensure the best performance of the equipment and to avoid dangerous situations:

- In daily use, choose the size that fits your fingers and avoid wearing it too loose or tight.
- Keep the heart rate blood oxygen sensor and lamp clean to prevent dirt from affecting the test data.
- If symptoms such as redness and swelling occur during use, discontinue use immediately and consult a physician.

Safety

- Keep the charging cable and device in a dry environment during charging.
- Do not touch the charging cable with wet hands.
- Keep the magnetic charging port surface of the ring clean to avoid poor contact and inability to charge.
- Do not allow water or other liquids to invade the charging cable.
- Do not disconnect or modify the battery, insert foreign objects, immerse the ring in water or other liquids for a long time, to avoid battery leakage, overheating, fire or explosion.
- Do not drop, squeeze or puncture the ring battery. Avoid exposing the battery to external pressures that may cause internal short-circuiting or overheating.